

Analysis of β ray spectroscopy

PH3105

2024-11-15

Debayan Sarkar
22MS002

Diptanuj Sarkar
22MS038

Sabarno Saha
22MS037



Contents

I. Introduction	1
II. Theory	1
II.1. β ray spectrometer	1
II.2. Operating principle of the β spectrometer	3
II.3. Decay Scheme of $^{22}_{11}\text{Na}$ and $^{90}_{38}\text{Sr}$	3
II.3.A. β decay	3
II.3.B. Decay scheme for $^{22}_{11}\text{Na}$	4
II.3.C. Decay scheme for $^{90}_{38}\text{Sr}$	4
II.4. Kinetic Energy of β Particles	5
III. Results and Analysis	6
III.1. $^{22}_{11}\text{Na}$ Source	6
III.1.A. Table for $^{22}_{11}\text{Na}$ Source	6
III.1.B. Plot of Energy vs Count for $^{22}_{11}\text{Na}$ Source	7
III.1.C. Energy corresponding to the maximum count for $^{22}_{11}\text{Na}$ source	7
III.2. $^{90}_{38}\text{Sr}$ Source	7
III.2.A. Table for $^{90}_{38}\text{Sr}$ Source	7
III.2.B. Plot of Energy vs Count for $^{90}_{38}\text{Sr}$ Source	9
III.2.C. Energy corresponding to the maximum count for $^{90}_{38}\text{Sr}$ source	10
IV. Sources of Error	10
V. Conclusion	10

I. Introduction

In this experiment, we perform β ray spectroscopy using a GM Counter. We measure the energy of the β particles emitted by $^{22}_{11}\text{Na}$ and $^{90}_{38}\text{Sr}$. We use this with a Hall probe to measure the energy of the β particles.

II. Theory

We lay out, in brief, the theory behind the β ray spectrometer used, and the decay of the radioactive sources that produce the β rays in our interaction.

II.1. β ray spectrometer




Figure 1 :: β spectrometry setup

A magnetic β ray spectrometer is an instrument that is designed to analyze the energy distribution of β particles emitted from a radioactive source. This type of spectrometer utilizes the phenomena of deflection of charged radiation particle in a magnetic field to discern the energy of the particles that are encountered by a Geiger-Mueller counter - a radiation detector that is not sensitive to the energy information of the incident radiation particles.

The crucial components of the β particle spectrometer are -

1. Power Supply: This is utilized in constant current mode to drive current through a solenoidal coil.
2. Solenoidal coil: The current flowing through the coil creates an approximately uniform magnetic field, whose magnitude can be controlled by varying the current allowed to flow through the coil. The magnitude of the magnetic field can be calculated using the following equation -

$$B = \mu n I$$

Where, B = Magnitude of the magnetic field. μ = Permeability of the material inside the solenoid. n = Number of turns in the solenoid per unit length of the solenoid. I = Current flowing through the solenoid.

1. Hall probe/Gaussmeter: It is used to measure the magnitude of the magnetic field present in the region where the β particles are present - using the Hall Effect. The Hall effect is a fundamental phenomenon observed in conductors (and semiconductors) when they are subjected to a magnetic field perpendicular to the direction of an electric current in the material. The effect of this is the generation of a transverse electric potential across the material, which is known as the Hall voltage. This is created due to the charge carriers in the material experiencing a transverse Lorentz force on them due to the magnetic field. This deflects the charge carriers moving under the influence of the electric field to one side of the conductor. This accumulation of charges creates an electric field transverse to the original electric field that was present in the material, which manifests in the form of a measurable voltage difference. The Hall voltage can be calculated using the following expression -

$$V_H = \frac{IB}{qnd} \quad [1]$$

V_H = Hall voltage I = Current flowing through the conductor (in Ampere) B = Magnitude of the magnetic field (in Tesla) q = charge of the carriers (in Coulomb) n = Number density of the charge carriers in the material d = Thickness of the conducting material. If calibrated correctly, the Hall Voltage is proportional to the magnetic field strength, and can thus be used to measure the magnetic field that is generated by the solenoid.

1. Circular aperture with diaphragms: This is the piece on top of the solenoid, which holds the Hall Probe, the radioactive source, and the Geiger-Mueller counter. The diaphragms present internally prevent the travel of any stray β particles from the source to the GM counter, except the particles that are travelling along a trajectory of a fixed radius.
2. Geiger-Mueller Counter: This is a widely used instrument for detecting and counting the number of ionizing radiation particles that are incident on it, including β rays. It consists of a cylindrical metal tube, a low pressure inert gas filling the tube, a quenching gas mixed in with the inert gas, and a central wire. The wall of the tube serves as the cathode, and the central wire works as the anode. A high voltage is applied between them. When a particle of ionizing radiation (a β particle in this case) enters the cylindrical chamber, it collides with the gas atoms that are inside. This interaction causes the ionization of the gas atoms, creating a positive ion and a free secondary electron. These free electrons are then accelerated by the electric field created in the tube by the application of the high voltage, and they drift towards the central wire. They accelerate and gain enough energy to collide with and ionize other gas atoms along the way, creating an avalanche like process where more and more atoms are ionised and a great amplification of charge takes place. This phenomenon is called Townsend Avalanche. These electrons are then picked up by the central wire, and the associated electronics detects the current pulse.

Each such detected pulse is counted as the incidence of one charged particle on the GM counter. It is to be noted that unlike ion chambers or proportional counters, these avalanches in the GM counters progress until a saturation level is reached - regardless of the energy of the initial incident radiation that triggered the events. In this process, the energy information of the incident radiation is lost entirely.

II.2. Operating principle of the β spectrometer

When the radioactive material in the source decays, it emits β particles of specific energies. These particles then travel into the annular top of the apparatus, where it is exposed to a magnetic field that is perpendicular to its direction of propagation in the plane. The β particles, which initially follow straight line trajectories, experience a Lorentz force given by -

$$(F) = q((v) \times (B)) \quad [2]$$

Where, (F) = Lorentz force (v) = Magnitude of the velocity of the β particle (B) = The magnetic field This force causes the particles to now move in a circular trajectory as it like a centripetal force, and the radius of the trajectory is given by -

$$R = \frac{p}{eB} \quad [3]$$

Where, R = Radius of the trajectory p = Momentum of the β particle e = Charge of the particle B = Magnitude of the magnetic field Due to the presence of the diaphragms in the annulus, all β particles except the ones present on trajectory of a predetermined radius are stopped from approaching the window of the GM counter. The GM counter is then used to count the number of particles that are incident over a fixed time period. Since the radius is fixed by the geometry of the equipment, $p \propto B$. Thus, by varying the current through the coil and measuring the magnitude of the magnetic field using the Hall probe, we can determine the energy of the β particles that can reach the GM counter. We count the number of such particles arriving at the counter, and we use this information to form the spectra of the sources.

II.3. Decay Scheme of $^{22}_{11}\text{Na}$ and $^{90}_{38}\text{Sr}$

We will detail out the β decay process and then We detail the β ray decay scheme of the $^{22}_{11}\text{Na}$ and $^{90}_{38}\text{Sr}$, both of which we will use in the experiment.

II.3.A. β decay

The β decay is a type of radioactive decay in which a beta particle is emitted. There are three types of beta decay processes: beta-minus decay, beta-plus decay or electron capture.

1. β^- decay : In beta-minus decay, a neutron decays into a proton, an electron, and an electron antineutrino.



2. β^+ decay : In beta-plus decay, a proton decays into a neutron, a positron, and an electron neutrino.



3. Electron Capture : In electron capture, an electron is captured by a proton in the nucleus, transforming it into a neutron and emitting a neutrino.



II.3.B. Decay scheme for $^{22}_{11}\text{Na}$

$^{22}_{11}\text{Na}$ undergoes a beta decay (β_+ emission) to form $^{22}_{10}\text{Ne}$ with an intermediate $^{22}_{10}\text{Ne}$ excited state. The excited state of $^{22}_{10}\text{Ne}$ then decays to the ground state by emitting a gamma ray. This is one of the three decay processes that take place. The other two are the electron capture and the gamma decay to the stable $^{22}_{10}\text{Ne}$. The gamma decay to $^{22}_{10}\text{Ne}$ occurs with negligible probability so we can ignore that. Here we consider the beta decay and the electron capture.



The decay scheme for $^{22}_{11}\text{Na}$ is shown below.






Figure 2 :: Decay scheme for $^{22}_{11}\text{Na}$, Source: Internet

II.3.C. Decay scheme for $^{90}_{38}\text{Sr}$

$^{90}_{38}\text{Sr}$ decays to form $^{90}_{40}\text{Zr}$ by beta decay (β_- emission) with an intermediate daughter nucleus $^{90}_{39}\text{Y}$. There are other decay processes. $^{90}_{39}\text{Y}$ can decay into an intermediate excited state of $^{90}_{40}\text{Zr}^*$. The decay scheme for $^{90}_{38}\text{Sr}$ is shown below.




and then $^{90}_{39}\text{Y}$ decays to $^{90}_{40}\text{Zr}^*$ by beta decay or



or directly having a β decay to $^{90}_{40}\text{Zr}$.



Figure 3 :: Decay scheme for $^{90}_{38}\text{Sr}$, Source: Internet

II.4. Kinetic Energy of β Particles

We use a Geiger-Muller counter to measure the energy of the gamma-ray emitted by the sources. Now, GM counters are not designed to measure energy since, in the Geiger Muller region of the characteristics curve, all energy information is lost. The spectrometer schema is laid out above. We find the relation between the counts and the energy of the beta-ray. Let us now derive the energy of the beta particle. The relativistic energy equation is given by:

$$E = \sqrt{p^2 c^2 + m_0 c^4} \quad [11]$$

where p is the momentum of the beta particle, m_0 is the rest mass of the beta particle, E is the total energy of the beta particle and c is the speed of light.

We also have $E = T + V = T + m_0 c^2$. Using this, we can write the above equation as:

$$T = \sqrt{p^2 c^2 + m_0 c^4} - m_0 c^2 \quad [12]$$

We place the beta source in a magnetic field and use the value for the average radius of curvature r to equate the Lorentz force with the centripetal force (mv^2/r). This gives us:

$$\frac{mv^2}{r} = q(E + vB) = qvB \Rightarrow p = mv = qBr \quad [13]$$

where q is the charge of the beta particle. Using this in the above equation, we get:

$$T = \sqrt{(qBr)^2 + m_0 c^4} - m_0 c^2 \quad [14]$$

III. Results and Analysis

III.1. $^{22}_{11}\text{Na}$ Source

For the $^{22}_{11}\text{Na}$ source, for each voltage we took 2 runs, for 60 seconds each. The energy was calculated using the formula derived above, in kilo electron volts.

III.1.A. Table for $^{22}_{11}\text{Na}$ Source

Magnetic Field(mT)	Voltage (V)	Current (A)	Counts (Run 1)	Counts (Run 2)	Avg. Counts	Energy (keV)
12.7	0.13	0.3	59	63	61	3.43E+01
25.9	0.26	0.6	116	102	109	1.31E+02
36.9	0.38	0.9	121	122	121	2.42E+02
39.4	0.41	1	125	121	123	2.70E+02
42.9	0.45	1.1	138	130	134	3.11E+02
47.7	0.5	1.2	113	130	121	3.68E+02
51.3	0.54	1.3	114	123	118	4.12E+02
55.6	0.58	1.4	108	104	106	4.67E+02
61.1	0.64	1.5	117	93	105	5.38E+02
69.5	0.74	1.8	66	76	71	6.50E+02
91.9	0.9	2.1	37	40	38	9.59E+02
95.8	0.99	2.4	49	42	45	1.01E+03

III.1.B. Plot of Energy vs Count for $^{22}_{11}\text{Na}$ Source




Figure 4 :: β -ray spectrum of $^{22}_{11}\text{Na}$

III.1.C. Energy corresponding to the maximum count for $^{22}_{11}\text{Na}$ source

The peak of the spectrum was observed to be at an energy of 310.36 keV with the average count of 134.

III.2. $^{90}_{38}\text{Sr}$ Source

For the $^{90}_{38}\text{Sr}$ source, we decided to take only one run for each voltage, again, for 60 seconds. However, this time the voltage was varied at a slower rate compared to that of the $^{22}_{11}\text{Na}$ source, to make sure we have more close-by datapoints to get a better and more accurate spectrum.

III.2.A. Table for $^{90}_{38}\text{Sr}$ Source

Magnetic Field(mT)	Voltage (V)	Current (A)	Counts	Corrected Counts (no background)	Energy (keV)
0	0	0	109	0	0.00E+00
4.6	0.06	0.1	216	107	4.63E+00
5.6	0.08	0.2	230	121	6.85E+00
7	0.1	0.2	256	147	1.07E+01
7.9	0.12	0.3	280	171	1.36E+01
9.4	0.14	0.3	301	192	1.91E+01
10.8	0.16	0.4	329	220	2.50E+01

12.2	0.18	0.4	327	218	3.18E+01
14.2	0.2	0.4	351	242	4.26E+01
15.1	0.22	0.5	398	289	4.79E+01
16.7	0.24	0.5	444	335	5.81E+01
17.9	0.26	0.6	428	319	6.62E+01
19	0.28	0.6	471	362	7.41E+01
20.6	0.3	0.7	488	379	8.61E+01
22.4	0.32	0.7	563	454	1.01E+02
24.6	0.34	0.8	609	500	1.19E+02
25.2	0.36	0.8	643	534	1.25E+02
27.3	0.38	0.9	666	557	1.44E+02
28.6	0.4	0.9	698	589	1.56E+02
30.7	0.42	1	772	663	1.77E+02
31.9	0.44	1	796	687	1.89E+02
33.6	0.46	1.1	826	717	2.07E+02
34.6	0.48	1.1	830	721	2.17E+02
36.7	0.5	1.2	951	842	2.40E+02
38.5	0.52	1.2	1034	925	2.60E+02
39.6	0.54	1.3	975	866	2.72E+02
41.9	0.56	1.3	1126	1017	2.99E+02
43.3	0.58	1.3	1126	1017	3.15E+02
45.4	0.6	1.4	1118	1009	3.40E+02
46.7	0.62	1.4	1256	1147	3.56E+02
48.4	0.64	1.5	1235	1126	3.77E+02
51.2	0.68	1.6	1353	1244	4.11E+02
54.7	0.72	1.7	1389	1280	4.55E+02
59	0.76	1.8	1480	1371	5.11E+02
62.2	0.8	1.9	1500	1391	5.52E+02
65.3	0.84	1.9	1545	1436	5.93E+02
70.5	0.9	2.1	1580	1471	6.63E+02
75.6	0.96	2.2	1547	1438	7.32E+02
81.2	1.02	2.4	1528	1419	8.09E+02
87.8	1.08	2.5	1494	1385	9.01E+02
92.2	1.14	2.7	1530	1421	9.63E+02
98.1	1.2	2.8	1414	1305	1.05E+03
102.2	1.26	2.9	1373	1264	1.10E+03
107.7	1.32	3.1	1265	1156	1.18E+03
113.3	1.38	3.3	1138	1029	1.26E+03
118.1	1.45	3.4	1081	972	1.33E+03

123.4	1.5	3.5	962	853	1.41E+03
128	1.56	3.7	869	760	1.47E+03
132.5	1.62	3.8	781	672	1.54E+03
137	1.68	4	748	639	1.61E+03
142.4	1.74	4.2	744	635	1.68E+03
149.1	1.8	4.3	561	452	1.78E+03
151.3	1.86	4.5	508	399	1.81E+03
155.9	1.92	4.6	508	399	1.88E+03
161.3	2	4.8	398	289	1.96E+03

III.2.B. Plot of Energy vs Count for $^{90}_{38}\text{Sr}$ Source




Figure 5 :: β -ray spectrum of $^{90}_{38}\text{Sr}$

III.2.C. Energy corresponding to the maximum count for $^{90}_{38}\text{Sr}$ source

The peak of the spectrum was observed to be at an energy of 662.83 keV with the average count of 1471.

IV. Sources of Error

1. A primary source of error is the comparatively high dead time of a Geiger-Mueller counter. This is a greater source of error for the Strontium source, due to the higher count rate. Any β radiation that is incident on the GM counter during the dead time will not trigger a series of Townsend avalanches and will not be registered as a count.
2. The efficiency of response of the GM counter is also a function of the energy of the incident radiation, and this introduces a systematic error into the spectrum that is being recorded.
3. Non-uniformity in the magnetic field created by passing current through the conducting coil can cause errors. Since the region where the β radiation is being emitted and captured lie outside of the solenoid, the field is not uniform. We have not accounted for the edge corrections in the magnetic field near the ends of a finite length solenoid, and thus these errors are present in the spectra that we have recorded.
4. Electrical fluctuations that might affect the counting equipment and the power sources used in the experiment are a source of error.
5. The β radiation particles also undergo some amount of energy straggling as they travel from the source to the GM counter, as the space between them is not an evacuated vacuum. It is reasonable to expect different β particles to have undergone different histories and thus have lost differing amounts of energy in the space between the detector and the source. This causes the spectra readings to be more diffuse than they ideally should be.
6. The assumption that the entire conical solid angle in which the radiation source puts out β particles has the same number density of the particles is erroneous, as we have not accounted for geometric factors such as the positioning of the radioactive material in the source. If the number density of β particles at some angle from the center of the cone (measured over a long enough period of time) changes as a function off the angle, this would introduce errors into the spectra that we have measured.

V. Conclusion

In this experiment we successfully recorded the β -ray spectra for $^{22}_{11}\text{Na}$ and $^{90}_{38}\text{Sr}$, using a β -ray spectrometer, and obtained the peak emission energy for both.