

LS2102: Biology Laboratory

LAB REPORT 2

Sabarno Saha
Roll No. 22MS037

April 25, 2024

Contents

1	Aim	2
2	Procedure	2
2.1	Preparation of Standard Curve	2
2.2	Estimation of Kinetic Parameters	2
3	Data Tables and Graphs	2
3.1	Data Table to Prepare Standard Absorbance vs pNP Curve	2
3.2	Table 2(Kinetic Parameters)	3
3.3	Table 3(Kinetic Parameters)	3
3.4	Graphs	4
4	Calculations and Results	6
5	Results	7

1 Aim

To determine the unknown concentration of p-NPP(para-Nitrophenylphosphate) and determine the enzymatic properties of Alkaline Phosphatase enzyme.

2 Procedure

2.1 Preparation of Standard Curve

1. Solutions of different concentrations (containing buffer, pNP and NaOH), having different amount of pNP were prepared.
2. Absorbance of each solution was taken at 405nm.
3. Graph of Absorbance vs. Concentration of pNP was plotted and fit with a $Ax + B$.
4. The value of slope and intercept from the fit curve was noted.

2.2 Estimation of Kinetic Parameters

1. The buffer, enzyme and substrate solutions were prepared according to the table provided.
2. The solutions were kept for incubation at 300K for 15 minutes.
3. Absorbance at 405nm was taken for each solution.
4. The reaction velocities were calculated by obtaining product at the specified absorbance from the standard linear fitting $Ax + b$ done in the previous step.
5. The Michaelis-Menton Curve was plotted with fit $\frac{Ax}{B + x}$
6. The Lineweaver-Burke Curve were plotted with again a standard linear fit
7. From the curve, various kinetic parameters were calculated.

3 Data Tables and Graphs

3.1 Data Table to Prepare Standard Absorbance vs pNP Curve

Tube Number	Volume of Buffer (μL)	Volume of pNP (μL)	Volume of NaOH (mL)	Concentration of pNP (μM)	Absorbance At 405 nm	Average Absorbance
Blank	1000	0	1	0	0.0	0.0
Blank					0.0	
1	985	15	1	3	0.024	0.021
2					0.018	
3	975	25	1	5	0.047	0.044
4					0.041	
5	950	50	1	10	0.075	0.0805
6					0.086	
7	925	75	1	15	0.131	0.147
8					0.163	
9	900	100	1	20	0.224	0.225
10					0.226	
11	875	125	1	25	0.272	0.273
12					0.274	
13	850	150	1	30	0.317	0.32
14					0.323	
15	800	200	1	40	0.41	0.417
16					0.424	


3.2 Table 2(Kinetic Parameters)

Tube No.	Volume of Buffer (μL)	Volume of pNP (μL)	Volume of ALP Enzyme (μL)	Volume of NaOH (mL)	Concentration of pNP (mM)	Absorbance At 405nm	Correction Factor	Corrected Absorbance
Blank	1000	0	25	1	0	0	0	0
Blank						0		0
1	965	10	25	1	0.1	0.021	0.001	0.02
2						0.012		0.011
3	950	25	25	1	0.25	0.038	0.004	0.034
4						0.037		0.033
5	925	50	25	1	0.5	0.057	0.010	0.047
6						0.061		0.051
7	900	75	25	1	0.75	0.056	0.011	0.045
8						0.060		0.049
9	875	100	25	1	1	0.076	0.014	0.062
10						0.080		0.066
11	775	200	25	1	2	0.117	0.029	0.088
12						0.124		0.095
13	725	250	25	1	2.5	0.142	0.040	0.102
14						0.138		0.098
15	475	500	25	1	5	0.209	0.074	0.135
16						0.215		0.141
17	225	750	25	1	7.5	0.278	0.112	0.2668
18						0.273		0.2618
19	75	900	25	1	9	0.326	0.133	0.193
20						0.334		0.201

3.3 Table 3(Kinetic Parameters)

Tube No.	Concentration of pNP (mM) [S]	1/[S] (mM ⁻¹)	Corrected Absorbance	Product from Standard Curve (μM)	V=P/t (μM/min)	Average (μM/min)	1/V (μM/min) ⁻¹	Average (μM/min) ⁻¹
1	0.1	10	0.02	2.684	0.178933	0.151367	5.588674	6.833109
2			0.011	1.857	0.1238		8.077544	
3	0.25	4	0.034	3.971	0.264733	0.261667	3.777386	3.822181
4			0.033	3.879	0.2586		3.866976	
5	0.5	2	0.047	5.165	0.344333	0.3566	2.904163	2.807585
6			0.051	5.533	0.368867		2.711007	
7	0.75	1.33	0.045	4.982	0.332133	0.344367	3.010839	2.907551
8			0.049	5.349	0.3566		2.804262	
9	1	1	0.062	6.544	0.436267	0.448533	2.292176	2.231157
10			0.066	6.912	0.4608		2.170139	
11	2	0.5	0.088	8.934	0.5956	0.617033	1.678979	1.622616
12			0.095	9.577	0.638467		1.566252	
13	2.5	0.4	0.102	10.22	0.681333	0.6691	1.46771	1.495045
14			0.098	9.853	0.656867		1.522379	
15	5	0.2	0.135	13.253	0.883533	0.901933	1.131819	1.109191
16			0.141	13.805	0.920333		1.086563	
17	7.5	0.13	0.2668	16.367	1.091133	1.109133	0.06109	0.060115
18			0.2618	16.907	1.127133		0.05914	
19	9	0.11	0.193	18.584	1.238933	1.263433	0.807146	0.791792
20			0.201	19.319	1.287933		0.776438	

3.4 Graphs




From fitting the line of the curve as per the equation $Ax+B$ we get the values of A,B as:-

B (y-intercept) = $-9.2011952191235e-03$

A (slope) = $1.0880478087649e-02$

Michelis-Menten Curve




From fitting the line of the curve as per the equation $\frac{Ax}{B+x}$ we get the values of A,B as:-

$$B = 2.45$$

$$A = 1.47$$

Linweaver-Burke Plot



From fitting the line of the curve as per the equation $Ax+B$ we get the values of A,B as:-

$$B \text{ (y-intercept)} = 1.19$$

$$A \text{ (slope)} = 0.597$$

4 Calculations and Results

From the graph we found out that

$$V_{max} = \frac{1}{B} = 0.840 \mu\text{M}/\text{min}.$$

$$K_{max} = A * V_{max} = 0.502 \text{ mM}$$

$$[E]_T = \frac{\text{Volume of Enzyme taken}}{\text{Total Volume}} \times \text{Concentration of Stock} = 0.504 \mu\text{M}$$

$$k_{cat} = 1.67 \text{ min}^{-1}$$

$$\text{Specific Activity of the enzyme} = \frac{0.840 * 1}{0.00025} = 3,360 \mu \text{ mol min}^{-1} \text{ mg}^{-1}$$

5 Results

Thus with the help of Michelis-Menten Curve and Linweaver-Burke Plot we could figure out the various kinetic parameters of the enzyme alkaline phosphate and the results are given below

$$V_{max} = 0.840\mu\text{M}/\text{min.}$$

$$K_{max} = 0.502 \text{ mM}$$

$$[E]_T = 0.504\mu\text{M}$$

$$k_{cat} = 1.67 \text{ min}^{-1}$$

$$\text{Specific Activity of the enzyme} = 3,360 \mu \text{ mol min}^{-1} \text{ mg}^{-1}$$